

Transformation of Yeast

1. Grow cells in liquid culture to mid-logarithmic stage (0.5-0.8 ODU/ml)
*cells at stationary phase have a thicker cell wall and are more difficult to transform
2. Centrifuge 5-10 ODU of cells at 5000 x *g* for 5 min at room temperature (RT)
3. Resuspend the cell pellet in 1 mL sterile water
4. Centrifuge at 5000 x *g* for 5 min at RT
5. Resuspend the cell pellet in 1 mL 0.1 M LiOAc/TE:
0.1 M lithium acetate
10 mM Tris pH 7.5
1 mM EDTA pH 8
6. Incubate in a 30°C water bath for 10 min
7. Centrifuge at 5000 x *g* for 5 min at RT
8. Resuspend the cell pellet in 50 µL 0.1 M LiOAc/TE
9. Add 10 µL of 10 mg/mL sheared salmon sperm DNA (stock is 10 mg/mL in sterile water at -20°C; thaw stock at 95°C for 5 min, chill on ice for 5 min, and vortex prior to using)
10. Add DNA to be transformed:
*if transforming closed circular miniprep DNA, use 1-5 µL
*if transforming linear PCR product for genomic integration, use the equivalent of 5-10 PCRs that have been precipitated and resuspended in 0.1 M LiOAc-TE
11. Add 700 µL of 40% PEG-4000 prepared in 0.1 M LiOAc/TE
*legend has it that the fresher the PEG solution, the higher the transformation efficiency
12. Incubate at 30°C for 30-60 min
13. Heat shock at 42°C for 20 min
14. Centrifuge at 5000 x *g* for 5 min at RT; aspirate supernatant
15. Resuspend cell pellet in 1 mL sterile water
16. Centrifuge at 5000 x *g* for 5 min at RT; aspirate supernatant
17. Resuspend cell pellet in 200 µL sterile water and spread solution onto YNB drop-out plate
18. Incubate until colonies appear (2-5 days depending upon the strain)

Note: The above protocol is relatively quick but does not yield as many colonies as is possible. To improve transformation efficiency (especially important for genomic integrations), incubate the cells overnight at 4°C during step 6 and overnight at 4°C again during step 12.

Also note: The above protocol is designed for transformation of strains that do not exhibit temperature-sensitive (ts) viability. If a ts strain is being transformed, adjust the incubation temperatures accordingly (you may want to eliminate the 42°C heat shock in step 13 if such a brief incubation at high temperature kills your strain).