Aridlands Ecology Lab Protocol Modified: 2009.08.05, S.Castle

Gravimetric Soil Moisture

Materials Needed:

Aluminum weigh-boats Sharpie Drying oven Field moist soil

Laboratory Processing:

- 1. Label the weigh-boat
- 2. Weigh empty boat and record boat weight
- 3. Tare and add ~10 g wet soil, record soil weight
- 4. Place samples in drying oven for 48 hours at 105°C
- 5. Take the final weight of dry soil plus tin

Calculation:

%Moisture = $\underbrace{(soil_{wet} - (soil_{drysoil+tin} - tin))}_{(soil_{drysoil+tin} - tin)*100}$