

Preparation of competent bacteria - Inoue method

Day 1:

1. Streak bacteria from frozen stock onto LB agar medium; incubate O/N at 37°C

Day 2:

2. Pick a single colony to inoculate 25 mL LB medium; shake at 37°C for 6-8 hr
3. Use starter culture to inoculate 3 250-mL LB medium (in 1-L flasks) *e.g., use 10 mL, 5 mL, and 3 mL (or 8 mL, 4 mL, and 2 mL)
4. Shake at 20°C O/N

Day 3:

5. Monitor OD₆₀₀ until it reaches 0.5, then chill the culture in ice water for 10 min
6. Spin out cells at 5,000 rpm for 10 min at 4°C
7. Decant and invert centrifuge tube onto paper towel for 2 min
8. Resuspend cell pellet in 80 mL ice-cold Inoue buffer
9. Repeat steps 6 and 7
10. Resuspend cell pellet in 20 mL ice-cold Inoue buffer and transfer volume to 50-mL disposable capped tube (e.g., 50-mL Falcon tube)
11. Add 1.5 mL DMSO; mix by inversion; incubate in ice water bath for 10 min
12. Dispense 210-µL aliquots into chilled eppi tubes (on ice); snap-freeze in liquid nitrogen.
13. Store at -80°C.

Typical competence is ~1x10⁶ colonies per µg DNA after transformation

Inoue Buffer	for 200 mL:	0.5M Pipes, pH 6.7
55 mM MnCl₂	2.17 g MnCl₂•4H₂O	15.1 g Pipes dissolved in 80 mL H₂O
15 mM CaCl₂	0.44 g CaCl₂•2H₂O	pH to 6.7 using 10M KOH
250 mM KCl	3.73 g KCl	Bring volume to 100 mL with H₂O
10 mM PIPES, pH 6.7	4 mL 0.5M PIPES, pH 6.7	