

Yeast Spore Enrichment for Random Spore Analysis

1. Start a 5ml culture of cells in Sporulation Media. Leave at RT in shaker for 1 week.
2. Prepare a suspension from $\sim 1 \times 10^8$ cells and asci of sporulated culture in a polypropylene 1.5ml μ fuge tube containing 180 μ l sterile DW. (1 ODU=3 X 10^7). Wash cells and asci 2X in 1ml sterile DW before suspending in 180 μ l DW.
3. Add 20 μ l of 5mg/ml solution of Zymolyase 20T in ZL buffer to the suspension and mix. Incubate at 30°C for 1 hour

| ZL Buffer | For 100ml |
|----------------------------------|-------------------------------------|
| 0.1 M NaPO ₄ , pH 6.5 | 10 ml 1M NaPO ₄ , pH 6.5 |
| 1.2 M Sorbital | 21.9g Sorbital |
| 40% Glycerol | 80 ml 50% Glycerol |
| DW | Up to 100ml |

Filter Steralize. (Allow time because thick)
Dissolve 5mg Zymolyase 20T in 1 ml ZL buffer
Store solution at -20°C

4. Spin 30 seconds FS, discard supernatant.
5. Suspend pellet in 1 ml DW and respin 30 seconds FS, discard supernatant.
6. Suspend pellet in 100 μ l DW and vortex on high speed for 2 minutes.
7. Pour the liquid out of the μ fuge tube and rinse the tube 2X with 1ml DW. Pour liquid out of tube each time. (Note: The hydrophobic spores will stick to the μ fuge tube wall while the cells and debris will remain mostly in suspension.)
8. Add 1.0ml sterile 0.01% NP-40 to the tube and sonicate for 30 seconds-1 minute.
9. Prepare 4X 1:10 serial dilutions of the sonicated suspension (50 μ l into 450 μ l 0.01% NP-40).
10. Spread 200 μ l of each dilution onto YPD plate or selective plate.