


# FUME HOOD VS



# BIOSAFETY CABINET



- Protects the user
- Used for chemical fumes/vapors
- Does not use a filter
- Exhausts air outside of the building



- Protects user, samples, & the lab environment
- Used for biologicals
- In most cases, filters air back into the room (no exhaust out of building)



## KNOW THE DIFFERENCE

**CU Green Labs Contact:**  
Kathy Ramirez-Aguilar, PhD  
[greenlabs@colorado.edu](mailto:greenlabs@colorado.edu)  
[colorado.edu/ecenter/greenlabs](http://colorado.edu/ecenter/greenlabs)  
303-735-5612



Green Labs Program  
UNIVERSITY OF COLORADO BOULDER  
In partnership with the Environmental Center, Environmental Health and Safety, and Facilities Management

