

# THINK!

## Before placing large lab equipment

???



Can the space handle the heat?

It's so hot in here!



**HELLO**  
my name is

**ULTRA LOW**

(so you can stop calling me minus 80)

Your samples may not require -80 °C freezing. Set your ultra low to the temperature your samples need, which may be -70 °C or higher.

## PLEASE: NO HALLWAYS OR SMALL ROOMS WITHOUT VENTILATION!

### Choose equipment rooms instead!

- Better for your equipment.
- Better for building ventilation energy savings.

CU Green Labs Contact:  
Kathy Ramirez-Aguilar  
[greenlabs@colorado.edu](mailto:greenlabs@colorado.edu)  
[ecenter.colorado.edu/greenlabs](http://ecenter.colorado.edu/greenlabs)  
303-492-8308



UNIVERSITY OF COLORADO BOULDER

Facilities Management  
Office of Sustainability  
Environmental Health and Safety  
Environmental Center

