

# ULT Freezer Emergency Response Plan



## Contact

<b>PI</b>		
Name	_____	
Phone	_____	
Email	_____	
<b>Lab Member</b>		
Name	_____	
Phone	_____	
Email	_____	
<b>Building Proctor</b>		
Name	_____	
Phone	_____	
Email	_____	



## Act

- Pack empty space in freezer with dry ice
- Contact CU Green Labs for a list of back up freezers near your lab
- Contact Ashlyn Norberg to borrow CU Green Labs mobile freezer

(ashlyn.norberg@colorado.edu or 301-802-9852)

Where	Hours	Contact
Facilities Management	8:00-4:30 Mon-Fri	Phone: (303) 492-6322 Fax: (303) 492-1186 Email: invorder@Colorado.EDU
MCDB Stores MCDB A3B40	8:00-4:30 Mon-Fri	Mike Retit Phone: (303) 492-6060
CHEM Stores Cristol Basement (Room #?)	8:00-4:30 Mon-Fri	Marshall Wilkenson Phone: (303) 492 6182
General Air 3750 Walnut St. Boulder CO	7:30-5:00 Mon-Fri	Phone: (303) 444-2138
King Soopers 1650 30 <sup>th</sup> St. Boulder CO	-24/7	Phone: (303) 443 9622



## Repair

1. Sienco Instrument Repair in Aurora, CO (Phil Bish, 303-934-1084, pjbish@centurylink.net)
2. Sercom Scientific Equipment Repair in Fort Collins and Wheat Ridge, CO (Ron Cummings, 303-356-1610, cummingsr@sercom-usa.com, Isabel Cruz, (303) 573-7266, cruzi@sercom-usa.com)

For advice on maintaining your ultra-low freezer or purchasing a new one, we invite you to contact Green Labs

CU Green Labs Contact:  
Kathy Ramirez  
greenlabs@colorado.edu  
ecenter.colorado.edu/greenlabs  
303-492-8308



UNIVERSITY OF COLORADO BOULDER

Facilities Management  
Office of Sustainability  
Environmental Health and Safety  
Environmental Center

