




RECYCLE LAB BROWN GLASS

Brown Glass ONLY

**NO
lids!!**



**250 mL
or
larger**

**No acid or base bottles
No oil, glycerol, or ethylene glycol bottles**

What To Do:

- ✓ Deface labels, write “Empty,” and PI’s name
- ✓ For Volatile Solvents:
Vent residual solvents from empty bottles in fume hoods for at least 24 hours
- ✓ For Solids and Non-Volatile Liquids:
Triple rinse empty bottles & dispose of rinse as hazardous waste
- ✓ Bottles should not have held any chemicals other than what is on the label
- ✓ Carefully place in totes to avoid breaking the glass

greenlabs@colorado.edu, 303-735-5612