

Available Rotors for the Shared Instrument Pool's Ultracentrifuges:

- *The rotors for the L8-70M centrifuges **CANNOT** be used with the Optima XE-90 centrifuge*
- *All of the older rotors have had their speed derated for safety and the maximum speed is written on each rotor*

Tabletop Ultracentrifuges		
Rotor Type	Number Available	Centrifuge
TLA 100.3	2	TL-100
TLA 120.2	1	TL-100
TLS 55	1	Optima MAX-TL

L8-70M Ultracentrifuges	
Rotor	Number Available
SW 41 Ti	1
SW 55 Ti	1
SW 28 Ti	1
Type 45 Ti	1
Type 70 Ti	2
Type 50.2 Ti	1

Optima XE-90 Ultracentrifuge	
Rotor	Number Available
SW 41 Ti	1
SW 32 Ti	1
SW 55 Ti	1
Type 45 Ti	1
Type 70 Ti	1
Type 50.2 Ti	1