

Available Rotors for the Shared Instrument Pool's High Speed Centrifuges:

- The rotors for the other centrifuges **CANNOT** be used with the Avanti JXN 26 centrifuges
- All of the older rotors have had their speed derated for safety and the maximum speed is written on each rotor

Beckman Avanti JXN 26 Super Speed Centrifuges	
Rotor	Number Available
JLA 12.500	2
JA-10	1
JLA-16.25	2
JA-25.50	2
JA-30.50 Ti	1

Beckman J2-21 Centrifuge	
Rotor	Number Available
JS-13.1	1
F10BA-6x500y	1
F14BA-6x250y	1
JA-20	1

Beckman Avanti J-20XP Centrifuge	
Rotor	Number Available
JLA-8.1	1
JA-20	1

Sorvall RC-5C Plus Centrifuges	
Rotor	Number Available
SS-34	2
SLA-3000	1
GS-3	1
GSA	1
SA-6000	2